

# Chemistry Lab 6 part II

Prerequisites:

- Mass of the aluminum used
- Petri dish with alum
- Percent yield of the product

## Melting point

1. Add some product to a capillary tube.
  - Gently tap the open end of the tube on your sample and push it in.
  - Tap the closed end of the tube to allow the sample to fall to the bottom of the tube
  - Fill 2–3 mm
2. Place the tube in the melting point apparatus. When you see melting, record the temperature. Once the entire sample is liquid, record the temperature.
  - These values give the "melting point range"
  - Tubes are waste

## Sulfate Ions

1. Dissolve a few crystals of alum in a test tube using 2–3 mL of water.
2. Add 2–3 drops of 1.0 M barium chloride
  - If there are sulfate ions, they will react with the barium

## Waters of Hydration

- We must determine how much mass is in the water sample
- Because our ionic compound has a far greater boiling point than water, we can heat the sample until the water boils away. (ensure that heat is the bare minimum to boil water and no greater)
- Given the initial mass and the final mass, we can then determine the amount of water present in the compound
- After a certain point, the mass of the sample will stop changing from heating, this is when the sample is dry

## Procedure help

- Ensure the crucible is clean and dry before use, then place over the Bunsen burner.
  - Lid should be ajar to allow steam to escape
- Only weight the crucible when cool, and heat it over a "cool" and low flame
- Turn of the burner when not in use
- Use ~1 gram of alum
- Stop heating immediately if white smoke or an odd smell is observed
- Use multiple heatings, and when the mass stops changing, the mission is complete

## Waste

alum salt is not dangerous, feel free to sink it.

Rinse the crucible before putting it away